



Infrared Blood Circulation Massager



Model: AB-753

Please read these instructions before use and keep for future reference

CONTENTS

Safety Instructions	1
Parts Identification and Control Panel	2
Remote	3
Operating Instructions	4
Maintenance & Storage	5
Specifications	5

SAFETY INSTRUCTIONS

When using an electrical appliance, basic precautions should always be followed, including the following:


- An appliance should never be left unattended when plugged in. Unplug from outlet when not in use and before putting on or taking off parts.
- Do not operate under flammable or flammable vapors. Excessive heating can occur and cause fire, electric shock, or injury to persons.
- Appliance operation is necessary when this appliance is used to, on or near clothes, linens, towels, bedding, pillows, or upholstered furniture.
- Use this appliance only for its intended use as described in the manual. Do not use attachments not recommended by the manufacturer.
- Never use the cord as a support for the appliance or plug. If it is not working properly, if it has been damaged or damaged, or dropped, stop using it. Contact your supplier for replacement or repair.
- Do not carry or plug the appliance into the power cord.
- Keep the cord away from heat, flames, and sharp edges.
- Do not use outdoors. This cord is not for household use only.
- To disconnect, turn the controls to the off position, then remove plug from the outlet.
- Do not use within four foot radius of other heat.
- Keep each session under 30 minutes.
- Do not hang on the cord.
- THIS UNIT IS NOT A MEDICAL DEVICE.

If you have the following conditions, please consult with your physician before using the unit:

1. Pregnant women or individuals who have undergone surgery and have inflamed or unhealed wounds.
2. Person who currently has an acute disease, has a history of heart problems, or has bone fractures that are not healed completely, or those with open sores or disease of the feet
3. Patient with loss of sensitivity to their lower limbs.
4. While using this unit, if you start to feel dizzy, lightheaded or uncomfortable, please stop the unit at once and consult your physician. Do not continue use until you have the permission of your doctor.

LEAVE THESE INSTRUCTIONS

PARTS IDENTIFICATION & CONTROL PANEL




1. Control panel
2. Message platform
3. CALORIE RESET button
4. POWER indicator
5. INFRARED button
6. SPEED buttons



7. TIMER indicators
8. AUTO mode indicator
9. POWER button
10. MANUAL mode indicator
11. TIMER button
12. CALORIE display

2

HEIMATTE



OPERATING INSTRUCTIONS

In traditional Chinese medicine practice, the feet are viewed as a mirror of your health. Massaging and stimulating the acupressure points will adjust the function of the entrails and bring your body back to its natural balance. Heated (infrared) treatment detoxifies and increases body's immune system. Unit is also suitable for use on various parts of the body.



1. Connect power cord to power source.
2. Place feet on the massaging nodes or place unit where desired.
3. Press POWER button to start. Unit begins in Auto mode. To select Manual mode, press POWER button on control panel.
4. To turn off, press Power button again.

AUTO FUNCTION

- The unit automatically switches between different levels of intensity and infrared treatment at 35 second intervals.
Low Intensity and speed level → increased level → increased level → increased level → decreased level → increased level
→ decreased level → decreased level → increased level → decreased level → increased level → increased level
- In Auto, intensity and infrared levels cannot be changed. The INFRARED and SPEED buttons are not valid
- To change massage time (5, 10 or 15 min), press the TIMER button on the control panel

MANUAL Function

- In Manual, you can choose from 8 infrared treatments and 12 vibrating speeds.
- Press the **INFRARED** button (on remote or control panel) to change treatment. LCD display will indicate level selected.
- Press the speed **▲** or **▼** button (on remote or control panel) to change speed. LCD display will indicate level selected.

4



Your Customer

If this product is found to be faulty as a result of faulty materials or workmanship within one year from date of purchase, it will be repaired free of charge.

This guarantee is subject to the following terms:

- Surpiston™ must be notified of the fault.
- Proof of purchase must be presented to Surpiston's nominated representative.
- The warranty will be void if the product is modified, misused or repaired by an unauthorised person.
- The warranty after repair will not be extended beyond the original one-year period.
- All replacement parts will be new or reconditioned.
- Parts, which are replaced, become the property of Surpiston.
- The warranty applies for the use of the product in the USA only.

What is NOT covered

- What is NOT COVERED:**
- Warranty does not include freight charges.
 - Incidental or consequential damage caused by possible defects with this product.
 - Damage to product caused by improper power supply voltage, accident, fire, floods or acts of nature.
 - Failure of product resulting from unauthorized modifications to the product.
 - Improper installation or failure to perform the necessary maintenance.

The GUARANTEE is in addition to your Statutory Rights.