

SPT®

Infrared Blood Circulation Massager



Model: AB-753

Please read these instructions before use and keep for future reference.

CONTENTS

Safety Instructions	1
Parts Identification and Control Panel	2
Remote	3
Operating Instructions	4-5
Maintenance & Storage	5
Specifications	5

SAFETY INSTRUCTIONS

When using an electrical appliance, basic precautions should always be followed, including the following:

1. An appliance should never be left unattended when plugged in. Unplug from outlet when not in use and before putting on or taking off parts.
2. Do not operate under blanket or pillow. Extreme heating can occur and cause fire, smoke attack, or injury to persons.
3. Close supervision is necessary when this appliance is used by, on, or near children, invalids, or disabled persons.
4. Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer. Never operate this appliance if it has a damaged cord or plug. If it is not working properly, it has been dropped or damaged, or dropped into water. Contact Surgeon-General for examination and repair.
5. Do not carry or drag the appliance by the power cord.
6. Keep the cord away from heated surfaces.
7. Do not use outdoors. This unit is for household use only.
8. To disconnect, turn all controls to the off position, then remove plug from the outlet.
9. Do not use within one hour before or after a meal.
10. Wipe each session under 30 minutes.
11. Do not stand on this device.
12. THIS UNIT IS NOT A MEDICAL DEVICE.

If you have the following conditions, please consult with your physician before using the unit:

1. Pregnant women or individuals who have undergone surgery and have reflexes or unsteady movements.
2. Person who currently has an acute disease, has a history of heart problems, or has slow heartbeat that is not under medical supervision, or those with open sores or disease of the feet.
3. Patient with loss of sensitivity to foot lower limbs.
4. While using this unit, if you start to feel dizzy, lightheaded or uncomfortable, please stop the unit at once and consult your physician. Do not continue use until you have the permission of your doctor.

SAVE THESE INSTRUCTIONS

1

PARTS IDENTIFICATION & CONTROL PANEL




1. Control panel
2. Massage platform
3. CALORE RESET button
4. POWER button
5. INFRARED button
6. SPEED buttons



7. TIMER indicators
8. AUTO mode indicator
9. POWER button
10. MANUAL mode indicator
11. TIMER button
12. CALORE buttons

2

REMOTE



3

OPERATING INSTRUCTIONS

In traditional Chinese medicine practice, the feet are viewed as a mirror of your health. Massaging and stimulating the acupoints on the feet will adjust the function of the organs and bring your body back to its natural balance. Treated (infrared) treatment decreases and increases body's immune system, and it also suitable for use on various parts of the body.



1. Connect power cord to power source.
2. Place feet on the massage platform or place unit where desired.
3. Press POWER button to start. Unit begins in Auto mode. To select Manual mode, press POWER button on control panel.
4. To turn off, press Power button again.

AUTO Function

- Unit automatically switches between different levels of intensity and infrared treatment at 30 second intervals.
- Low infrared and speed level → increased level → increased level → increased level → decreased level → decreased level → decreased level → increased level → increased level → increased level → increased level → increased level → increased level.
- In Auto, intensity and infrared levels cannot be changed. The INFRARED and SPEED buttons are not valid.
- To change massage time (5, 10 or 15 min), press the TIMER button on the control panel.

MANUAL Function

- In Manual, you can choose from 8 infrared treatments and 12 vibrating speeds.
- Press the INFRARED button (on remote or control panel) to change treatment. LCD display will indicate level selected.
- Press the speed ▲ or ▼ button (on remote or control panel) to change speed. LCD display will indicate level selected.

4

- Press TIMER button to set massage time (5, 10 or 15 minutes).
- CALORE count and display. Once unit begins operating, the LCD display will include a "Calore Count". The number displayed is 17100 and is for general reference only. It does not represent actual calories burned. To read Calore count display, press the CALORE RESET button on the control panel.

MAINTENANCE & STORAGE

- Turn power off and unplug from power source before maintenance.
- Do not spray water or moisture onto or any parts of it in water or any other type of liquid.
- To clean, wipe with damp cloth.
- For storage, it is suggested to place unit back in its original packaging and store in a cool, dry place.

SPECIFICATIONS

Rated voltage	~120V / 60Hz
Rated power	60W
Cord dimension (L x W x H)	145.5 x 12 x 2.75 in
Package dimension (L x W x H)	16.5 x 14 x 9 in
Net weight (main body only)	5.6 lbs
Gross weight	7.5 lbs

5

SPT
Your Satisfaction

If this product is found to be faulty as a result of faulty materials or workmanship within one year from date of purchase, it will be repaired free of charge.

This guarantee is subject to the following terms:

- Replacement must be called at the factory.
- Proof of purchase must be presented to Surgeon-General's nominated representative.
- The warranty will be void if the product is modified, misused or repaired by an unauthorized person.
- The warranty after repair will not be extended beyond the original one-year period.
- All replacement parts will be new or reconditioned.
- Parts, which are replaced, become the property of Surgeon-General.
- The warranty applies for the use of the product in the USA only.

What is NOT COVERED:

- Warranty does not include freight charges.
- Incidental or consequential damage caused by possible defects with this product.
- Damage to product caused by improper power supply voltage, accident, fire, floods or acts of nature.
- Failure of product resulting from unauthorized modifications to the product.
- Improper installation or failure to perform the necessary maintenance.

This GUARANTEE is in addition to your Statutory Rights